

CLO-TABS

Convenient and easy to use economical chlorine tablets



CLO-TABS are convenient chlorine release tablets designed for simple dilution with water to provide a fresh chlorine solution of known strength.

ADVANTAGES

Convenient stable chlorine release tablets in a re-sealable tub.

Provides accurate dilution and reduces waste.

1 tablet dissolved in 5 litres water produces a 200ppm available chlorine solution.

Economical in use.

Long life, tablets maintain their strength during storage.

Reduced storage space and easy to handle.

Exceptional stain and tannin removal.

NSF approved

Powerful disinfectant.

TYPICAL APPLICATIONS

- Use as a disinfectant for toilets urinals and sinks.
- Disinfects and deodorises mop heads.
- Bleaches and removes stains from dishcloths and tea towels
- Ideal for removing tannin and stains from cutlery, cups and chopping boards.
- Suitable for fruit, salad and vegetable washing.

APPLICATION INSTRUCTIONS

1. Dissolve the appropriate number of tablets in water to give the strength and volume of solution required.
 2. Allow tablets to disperse and then mix the solution thoroughly prior to use.
 3. In order to ensure that chlorine levels are maintained it is recommended for critical applications that stock solutions are replenished daily.
- For general disinfection of food contact surfaces, equipment, walls, floors and equipment following cleaning. Apply Clo-Tab solution; allow 5 minutes contact time and rinse thoroughly with clean water.
 - For fruit salad and vegetable washing rinse of excess soil. Allow to soak for 5 minutes and then rinse thoroughly with fresh clean water.
 - When removing tannin stains allow articles to soak for up to 15 minutes and then rinse thoroughly with clean water.
 - For soak cleaning of cloths and mops allow to soak for 1-2 hours and then rinse through thoroughly with clean water.
 - To disinfect drains and waste pipes pour the Clo-Tab solution into the drain and allow to stand for up to 30 minutes. Flush or rinse with water.



NB. Do not mix Clo-Tab solutions with acids or any other chemical.

CLO-TABS

Convenient and easy to use economical chlorine tablets

Application	Recommended Chlorine content	Tablet dilution
General disinfection of food contact surfaces	200ppm	1 tablet in 5lts of water
Fruit salad and vegetable washing	Light soiling 75ppm Heavy soiling 125ppm	1 tablet in 13 lts of water 1 tablet in 8 lts of water
Tannin Removal	200-400ppm	1-2 tablets in 5 lts of water
Soak cleaning of cloths and mops	400ppm	2 tablets in 5 lts of water
Disinfection of drains	200-1000ppm	1-5 tablets in 5 lts of water
Disinfection of urine and faeces	1000ppm	5 tablets in 5 lts of water

TECHNICAL DATA

APPEARANCE	Tablet
COLOUR	White
ODOUR	Chlorine
SOLUBILITY	Soluble in water
pH VALUE (conc. solution)	6.0 @ 1.0%
PACK SIZES	➤ (6 x 200 3.25g Tablets)

HEALTH AND SAFETY, STORAGE AND PRECAUTIONS

Refer to Safety Data Sheet (available on request)

Use biocides safely. Always read the label and product information before use.

ACCREDITATIONS & REGISTRATIONS

➤ NSF (D2) 160742

ALSO AVAILABLE



BIOCIDAL CLEANER

Biocidal hard surface cleaner



HYPOSAN

High strength hypochlorite bleach for general purpose sanitising and destaining



WRS3 WASHROOM CLEANER CONCENTRATE

Concentrated daily washroom cleaner and descaler for use on all washroom surfaces



Manufactured by: Arrow Solutions, Rawdon Road, Moira, Swadlincote, Derbyshire, DE12 6DA, UK.
+44 (0)1283 221044 . sales@arrowchem.com . arrowchem.com

A division of Reabrook Ltd.
Registered in England No.804733